


Anti-Fen-1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - P39748

Item Number: CM0004600



Introduction

FEN-1 Rabbit Monoclonal Antibody validated for WB, IHC-P, and ELISA. Detects endogenous levels of FEN-1 in Human and Mouse. FEN-1 is a key nuclease in DNA replication and repair.

[Learn More](#)

Parameter	Value
Product Name	FEN-1 Rabbit mAb
Remarks/Alias	MF1; RAD2; FEN-1
Species	Human
Gene ID (Human)	2237
Gene ID	2237
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 250-350 of human FEN-1 (P39748).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P39748
Protein Weight	43kDa
Shipping	Ice bag