

## Anti-Fer Polyclonal Antibody For Wb And Elisa - P16591

Item Number: CM0030295



### Introduction

Rabbit polyclonal antibody targeting human FER (TYK3/p94-Fer) for Western blot and ELISA. Detects FER in human, mouse, and rat samples. Ideal for studying cell adhesion, migration, and receptor signaling.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | FER Rabbit pAb                                                                                                   |
| Remarks/Alias    | TYK3; PPP1R74; p94-Fer; FER                                                                                      |
| Species          | Human                                                                                                            |
| Gene ID (Human)  | 2241                                                                                                             |
| Gene ID          | 2241                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 200-450 of human FER (NP_005237.2) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                                |
| SWISS            | P16591                                                                                                           |
| Protein Weight   | 95kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |