

## Anti-Ferredoxin Reductase Polyclonal Antibody For Wb, If-P, Elisa - P22570

Item Number: CM0019837

### Introduction

Rabbit polyclonal antibody targeting Ferredoxin Reductase (FDR). Suitable for WB, IF-P, ELISA. Cross-reacts with human, mouse, rat. Investigate mitochondrial P450 electron transfer and coenzyme Q biosynthesis.

[Learn More](#)

| Parameter        | Value                                                                                                                             |
|------------------|-----------------------------------------------------------------------------------------------------------------------------------|
| Product Name     | Ferredoxin Reductase Rabbit pAb                                                                                                   |
| Remarks/Alias    | ADR; ADXR; ANOA; FDR                                                                                                              |
| Species          | Human                                                                                                                             |
| Gene ID (Human)  | 2232                                                                                                                              |
| Gene ID          | 2232                                                                                                                              |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 142-491 of human Ferredoxin Reductase (NP_077728.3) |
| Source           | Rabbit                                                                                                                            |
| Category         | Polyclonal Antibodies                                                                                                             |
| Application      | WB, IF-P, ELISA                                                                                                                   |
| Cross Reactivity | Human, Mouse, Rat                                                                                                                 |
| SWISS            | P22570                                                                                                                            |
| Protein Weight   | 54kDa                                                                                                                             |
| Shipping         | Ice bag                                                                                                                           |