

Anti-Ferredoxin Reductase Rabbit Monoclonal Antibody For Wb - P22570

Item Number: CM0009193

Introduction

Rabbit monoclonal antibody targeting human Ferredoxin Reductase (FDXR). Validated for WB, IF/ICC, IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Suitable for mitochondrial P450 system and steroid metabolism research. Predicted band 54kDa.

[Learn More](#)

Parameter	Value
Product Name	Ferredoxin Reductase Rabbit mAb
Remarks/Alias	ADR; ADXR; ANOA
Species	Human
Gene ID	2232
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 142-491 of human Ferredoxin Reductase (NP_077728.3)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22570
Protein Weight	54kDa
Shipping	Ice bag