


Anti-Ferritin Heavy Chain Monoclonal Antibody For Wb, If-P, Elisa - P02794

Item Number: CM0002721



Introduction

Recombinant rabbit monoclonal antibody targeting human Ferritin Heavy Chain (FTH1). Validated for Western blot, immunofluorescence, and ELISA; cross-reacts with mouse and rat. Ideal for iron metabolism, autophagy, and neurodegeneration research.

[Learn More](#)

Parameter	Value
Product Name	Ferritin Heavy Chain Rabbit mAb
Remarks/Alias	FHC; FTH; HFE5; PLIF; FTHL6; PIG15; Ferritin Heavy Chain
Species	Human
GeneID (Human)	2495
GeneID	2495
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Ferritin Heavy Chain (P02794).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02794
Protein Weight	21kDa
Shipping	Ice bag