

Anti-Ferritin Heavy Chain Monoclonal Antibody For Wb, If-P, Ihc-P, Elisa - P02794

Item Number: CM0004856

Introduction

Ferritin Heavy Chain rabbit monoclonal antibody validated for WB, IF-P, IHC-P, ELISA in human, mouse, rat. Targets FTH1, a key iron-storage protein with ferroxidase activity. Suitable for autophagy, iron homeostasis, and neurodegeneration studies.

[Learn More](#)

Parameter	Value
Product Name	Ferritin Heavy Chain Rabbit PolymAb®
Remarks/Alias	FHC; FTH; HFE5; PLIF; FTHL6; PIG15; Ferritin Heavy Chain
Species	Human
GeneID (Human)	2495
GeneID	2495
Immunogen	Synthetic peptide derived from human Ferritin Heavy Chain
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02794
Protein Weight	21kDa
Shipping	Ice bag