


Anti-Ferritin Light Chain Polyclonal Antibody For Wb, Ihc-P, Elisa - P02792

Item Number: CM0019911



Introduction

Rabbit polyclonal antibody targeting human Ferritin Light Chain (FTL). Validated for WB, IHC-P, ELISA. Cross-reacts with mouse and rat. Useful for iron homeostasis and neurodegeneration research.

[Learn More](#)

Specification	Details
Product Name	Ferritin Light Chain Rabbit pAb
Remarks/Alias	LFTD; NBIA3; Ferritin Light Chain
Species	Human
GeneID (human)	2512
GeneID	2512
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-175 of human Ferritin Light Chain (NP_000137.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P02792
Protein Weight	20 kDa
Shipping	Ice bag