


Anti-Ferrochelatase (Fech) Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - P22830

Item Number: CM0030742



Introduction

Rabbit polyclonal antibody targeting human ferrochelatase (FECH), suitable for Western blot, immunohistochemistry (IHC-P), and ELISA. Detects human, mouse, and rat FECH, supporting heme biosynthesis and mitochondrial studies.

[Learn More](#)

Parameter	Value
Product Name	FECH Rabbit pAb
Remarks/Alias	EPP; FCE; EPP1; FECH
Species	Human
Gene ID	2235
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human FECH (NP_000131.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22830
Protein Weight	48kDa
Shipping	Ice bag