

Anti-Fetuin A/Ahsg Monoclonal Antibody For Wb, Elisa - P02765

Item Number: CM0002107

Introduction

Rabbit monoclonal antibody targeting human Fetuin A/AHSG, validated for WB and ELISA applications. Ideal for studies of bone mineralization, endocytosis, and plasma protein function. SWISS: P02765.

[Learn More](#)

Parameter	Value
Product Name	Fetuin A/AHSG Rabbit mAb
Remarks/Alias	AHS; A2HS; HSGA; APMR1; FETUA; Fetuin A/AHSG
Species	Human
GeneID (Human)	197
GeneID	197
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Fetuin A/AHSG (P02765)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02765
Protein Weight	39kDa
Shipping	Ice bag