

Anti-Fetuin A/Ahsg Polyclonal Antibody For Wb, Elisa - P02765

Item Number: CM0019190

Introduction

Rabbit polyclonal antibody specific to human Fetuin A/AHSG for Western blotting and ELISA. Ideal for research in bone mineral homeostasis, endocytosis, and inflammation. - P02765

[Learn More](#)

Parameter	Value
Product Name	Fetuin A/AHSG Rabbit pAb
Remarks/Alias	AHS; A2HS; HSGA; APMR1; FETUA; Fetuin A/AHSG
Species	Human
Gene ID (Human)	197
Gene ID	197
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 19-300 of human Fetuin A/Fetuin A/AHSG (NP_001613.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P02765
Protein Weight	39kDa
Shipping	Ice bag