


Anti-Ffar3 Rabbit Polyclonal Antibody For Wb, Elisa - O14843

Item Number: CM0016334



Introduction

Rabbit polyclonal antibody against human FFAR3 (GPR41), validated for WB and ELISA, cross-reacts with mouse and rat. For research on energy homeostasis, gut immunity, and short-chain fatty acid signaling.

[Learn More](#)

Product Name	FFAR3 Rabbit pAb
Remarks/Alias	FFA3R; GPR41; GPR42; FFAR3
Species	Human
GeneID (Human)	2865
GeneID	2865
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 257-346 of human FFAR3 (NP_005295.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O14843
Protein Weight	39kDa
Shipping	Ice bag