

Anti-Ffar4 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Nul3

Item Number: CM0023620



Introduction

High-quality rabbit polyclonal antibody targeting human FFAR4 (GPR120), validated for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Ideal for research on fatty acid sensing, adipogenesis, and metabolic inflammation.

[Learn More](#)

Parameter	Value
Product Name	FFAR4 Rabbit pAb
Remarks/Alias	GT01; PGR4; OB10Q; BMIQ10; GPR120; GPR129; O3FAR1; FFAR4
Species	Human
GeneID (Human)	338557
GeneID	338557
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 100-200 of human FFAR4 (NP_859529.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5NUL3
Protein Weight	40kDa
Shipping	Ice bag