

Anti-Ffar4 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Nul3

Item Number: CM0020270



Introduction

FFAR4 Rabbit Polyclonal Antibody (pAb) detects human, mouse, rat FFAR4 (GPR120) in WB and ELISA. Ideal for metabolism, inflammation, and adipogenesis research.

[Learn More](#)

Parameter	Value
Product Name	FFAR4 Rabbit pAb
Remarks/Alias	GT01; PGR4; OB10Q; BMIQ10; GPR120; GPR129; O3FAR1; FFAR4
Species	Human
Gene ID (Human)	338557
Gene ID	338557
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 333-377 of human FFAR4 (NP_859529.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q5NUL3
Protein Weight	40kDa
Shipping	Ice bag