

Anti-Fgd1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P98174

Item Number: CM0022888

Introduction

Rabbit polyclonal antibody against human FGD1 (FYVE, RhoGEF and PH domain-containing protein 1), validated for WB, IHC-P, IF/ICC, and ELISA. Suitable for human, mouse, and rat samples. Ideal for studying CDC42-mediated actin cytoskeleton regulation and cell morphology.

[Learn More](#)

Parameter	Value
Product Name	FGD1 Rabbit pAb
Remarks/Alias	AAS; FGDY; MRXS16; ZFYVE3; FGD1
Species	Human
GeneID	2245
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 700-800 of human FGD1 (NP_004454.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P98174
Protein Weight	107kDa
Shipping	Ice bag