

Anti-Fgf1 Polyclonal Antibody For If/Icc, Elisa - P05230

Item Number: CM0016384



Introduction

Anti-FGF1 Rabbit Polyclonal Antibody validated for IF/ICC and ELISA. Cross-reacts with Human, Mouse, Rat. Investigates FGF1's role in cell survival, angiogenesis, and differentiation.

[Learn More](#)

Parameter	Value
Product Name	FGF1 Rabbit pAb
Remarks/Alias	AFGF; ECGF; FGFA; ECGFA; ECGFB; FGF-1; HBGF1; HBGF-1; GLIO703; ECGF-beta; FGF-alpha; FGF1
Species	Human
GeneID (Human)	2246
GeneID	2246
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05230
Protein Weight	17kDa
Shipping	Ice bag