

Anti-Fgf1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P05230

Item Number: CM0017319



Introduction

Rabbit polyclonal antibody against human FGF1, validated for WB, IHC-P, IF/ICC, ELISA. Cross-reacts with mouse and rat. Suitable for research into cell growth, angiogenesis, and differentiation.

[Learn More](#)

Parameter	Value
Product Name	FGF1 Rabbit pAb
Remarks/Alias	AFGF; ECGF; FGFA; ECGFA; ECGFB; FGF-1; HBGF1; HBGF-1; GLIO703; ECGF-beta; FGF-alpha; FGF1
Species	Human
Gene ID (Human)	2246
Gene ID	2246
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 16-155 of human FGF1 (NP_000791.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P05230
Protein Weight	17kDa
Shipping	Ice bag