

## Anti-Fgf10 Rabbit Polyclonal Antibody For Wb, Elisa - O15520

Item Number: CM0016036



### Introduction

Anti-FGF10 rabbit polyclonal antibody validated for Western blot (WB) and ELISA. Cross-reacts with mouse and rat. Ideal for studying embryonic development, cell proliferation, and differentiation.

[Learn More](#)

| Parameter        | Value                 |
|------------------|-----------------------|
| Product Name     | FGF10 Rabbit pAb      |
| Remarks/Alias    | LADD3; FGF10          |
| Species          | Human                 |
| GeneID (Human)   | 2255                  |
| GeneID           | 2255                  |
| Immunogen        | Recombinant protein   |
| Source           | Rabbit                |
| Category         | Polyclonal Antibodies |
| Application      | WB, ELISA             |
| Cross Reactivity | Mouse, Rat            |
| SWISS            | O15520                |
| Protein Weight   | 23kDa                 |
| Shipping         | Ice bag               |