


Anti-Fgf13 Rabbit Polyclonal Antibody For Wb - Q92913

Item Number: CM0025174



Introduction

Rabbit polyclonal antibody against human FGF13, validated for WB and ELISA. Detects human, mouse, rat FGF13, involved in neuronal development and sodium channel regulation. Ideal for neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	FGF13 Rabbit pAb
Remarks/Alias	FGF2; FHF2; DEE90; FHF-2; FGF-13; XLID110; LINC00889; FGF13
Species	Human
GeneID (Human)	2258
GeneID	2258
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-245 of human FGF13 (NP_004105.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q92913
Protein Weight	28kDa
Shipping	Ice bag