

## Anti-Fgf17 Polyclonal Antibody For Wb, If/Icc, Elisa - O60258

Item Number: CM0019460

### Introduction

Anti-FGF17 Rabbit Polyclonal Antibody validated for WB, IF/ICC, and ELISA. Detects human, mouse, and rat FGF17, a secreted growth factor crucial for embryonic brain development. Ideal for developmental biology studies.

[Learn More](#)

| Parameter        | Value                       |
|------------------|-----------------------------|
| Product Name     | FGF17 Rabbit pAb            |
| Remarks/Alias    | HH20; FGF-13; FGF-17; FGF17 |
| Species          | Human                       |
| GenelD (human)   | 8822                        |
| GenelD           | 8822                        |
| Immunogen        | Recombinant protein         |
| Source           | Rabbit                      |
| Category         | Polyclonal Antibodies       |
| Application      | WB, IF/ICC, ELISA           |
| Cross Reactivity | Human, Mouse, Rat           |
| SWISS            | O60258                      |
| Protein Weight   | 25kDa                       |
| Shipping         | Ice bag                     |