

Anti-Fgf2 Rabbit Monoclonal Antibody For Wb, Elisa - P09038

Item Number: CM0007034



Introduction

Rabbit monoclonal antibody targeting human FGF2, suitable for Western blot and ELISA. Detects multiple isoforms (17-30 kDa). Cross-reacts with mouse and rat. Useful for studying cell survival, angiogenesis, and differentiation.

[Learn More](#)

Parameter	Value
Product Name	FGF2 Rabbit mAb
Remarks/Alias	BFGF; FGFB; FGF-2; HBGF-2; FGF2
Species	Human
Gene ID (Human)	2247
Gene ID	2247
Immunogen	Synthetic peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P09038
Protein Weight	17kDa/21kDa/22kDa/30kDa
Shipping	Ice bag