


Anti-Fgf21 Polyclonal Antibody For Wb, Elisa - Q9Nsa1

Item Number: CM0017552



Introduction

Rabbit polyclonal antibody against human FGF21, validated for WB and ELISA. Detects human and mouse FGF21, a secreted growth factor regulating glucose homeostasis and insulin sensitivity.

[Learn More](#)

Parameter	Value
Product Name	FGF21 Rabbit pAb
Remarks/Alias	FGF21
Species	Human
GeneID (Human)	26291
GeneID	26291
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 29-209 of human FGF21 (NP_061986.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9NSA1
Protein Weight	22 kDa
Shipping	Ice bag