


Anti-Fgf21 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Nsa1

Item Number: CM0006417



Introduction

Rabbit monoclonal antibody targeting human fibroblast growth factor 21 (FGF21). Validated for WB, IHC-P, ELISA. Recognizes human, mouse, and rat FGF21. Ideal for metabolic and glucose homeostasis research.

[Learn More](#)

Parameter	Specification
Product Name	FGF21 Rabbit mAb
Remarks/Alias	FGF21; fibroblast growth factor 21
Species	Human
GeneID (Human)	26291
GeneID	26291
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 1-100 of human FGF21 (NP_061986.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NSA1
Protein Weight	22kDa
Shipping	Ice bag