

Anti-Fgf23 Rabbit Monoclonal Antibody For Wb, If/Icc, Elisa - Q9Gzv9

Item Number: CM0007624



Introduction

FGF23 rabbit monoclonal antibody validated for WB, IF/ICC, and ELISA. Detects human FGF23, a key regulator of phosphate and vitamin D metabolism. Ideal for bone and kidney disease research.

[Learn More](#)

Parameter	Value
Product Name	FGF23 Rabbit mAb
Remarks/Alias	ADHR; FGFN; HYPF; HFTC2; HPDR2; PHPTC; FGF23
Species	Human
GeneID (human)	8074
GeneID	8074
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-251 of human FGF23 (NP_065689.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human
SWISS	Q9GZV9
Protein Weight	28kDa
Shipping	Ice bag