

Anti-Fgf3 Polyclonal Antibody For Wb, Elisa - P11487

Item Number: CM0021070

Introduction

FGF3 Rabbit Polyclonal Antibody validated for WB and ELISA. Detects human FGF3 and cross-reacts with rat. Ideal for embryonic development and cell differentiation research.

[Learn More](#)

Parameter	Value
Product Name	FGF3 Rabbit pAb
Remarks/Alias	INT2; HBGF-3; FGF3
Species	Human
GeneID (Human)	2248
GeneID	2248
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P11487
Protein Weight	27kDa
Shipping	Ice bag