


Anti-Fgf4 Rabbit Polyclonal Antibody For Wb And Elisa - P08620

Item Number: CM0016148



Introduction

High-quality rabbit polyclonal antibody targeting human fibroblast growth factor 4 (FGF4).

Validated for Western blot and ELISA applications, with cross-reactivity to rat. Ideal for embryonic development, cell differentiation, and growth factor research.

[Learn More](#)

Parameter	Value
Product Name	FGF4 Rabbit pAb
Remarks/Alias	HST, KFGF, FGF-4, HST-1, HSTF1, K-FGF, HBGF-4, HSTF-1, FGF4
Species	Human
Gene ID	2249
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 31-206 of human FGF4 (NP_001998.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Rat
SWISS	P08620
Protein Weight	22kDa
Shipping	Ice bag