

Anti-Fgf6 Polyclonal Antibody For Wb, Elisa - P10767

Item Number: CM0030563

Introduction

Affinity-purified rabbit anti-FGF6 polyclonal antibody for WB, ELISA. Targets human FGF6 (HST2/HBGF-6); cross-reacts with mouse and rat. Ideal for cell proliferation, differentiation, angiogenesis, and myogenesis studies.

[Learn More](#)

Parameter	Value
Product Name	FGF6 Rabbit pAb
Remarks/Alias	HST2; HBGF-6; FGF6
Species	Human
GeneID (Human)	2251
GeneID	2251
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGF6 (NP_066276.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10767
Protein Weight	23kDa
Shipping	Ice bag