

Anti-Fgf7/Kgf Rabbit Monoclonal Antibody For Wb, Elisa - P21781

Item Number: CM0009057

Introduction

Recombinant rabbit monoclonal antibody against human FGF7/KGF, validated for Western blot and ELISA. Reacts with human, mouse, and rat. Suitable for research on growth factors and epithelial cell proliferation.

[Learn More](#)

Parameter	Details
Product Name	FGF7/KGF Rabbit mAb
Remarks/Alias	KGF; HBGF-7; Fibroblast growth factor 7; FGF7/KGF
Species	Human
Gene ID (Human)	2252
Gene ID	2252
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGF7 (NP_002000.1)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21781
Protein Weight	12kDa/23kDa
Shipping	Ice bag