


Anti-Fgf9 Polyclonal Antibody For Wb And Elisa - P31371

Item Number: CM0028679



Introduction

Rabbit polyclonal antibody targeting FGF9 (GAF, FGF-9, SYNS3, HBFG-9) for Western blot and ELISA. Validated in human, mouse, rat. Ideal for developmental biology, glial cell, and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	FGF9 Rabbit pAb
Remarks/Alias	GAF; FGF-9; SYNS3; HBFG-9; HBGF-9; FGF9
Species	Human
GeneID (human)	2254
GeneID	2254
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P31371
Protein Weight	23kDa
Shipping	Ice bag