

Anti-Fgfr1Op Polyclonal Antibody For Wb, Elisa - O95684

Item Number: CM0023727



Introduction

High-quality FGFR1OP rabbit polyclonal antibody validated for Western blot and ELISA. Reacts with human, mouse, and rat samples. Targets centrosomal protein 43 for cytoskeleton and cilia research.

[Learn More](#)

Parameter	Value
Product Name	FGFR1OP Rabbit pAb
Remarks/Alias	FOP; FGFR1OP
Species	Human
GeneID (human)	11116
GeneID	11116
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-70 of human FGFR1OP (NP_919410.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O95684
Protein Weight	43kDa
Shipping	Ice bag