

Anti-Fgfr2 (Ko Validated) Monoclonal Antibody For Wb, Elisa - P21802

Item Number: CM0002436



Introduction

Rabbit monoclonal antibody targeting human FGFR2 (P21802), validated for Western blot (WB) and ELISA. Suitable for research in cell signaling, development, and cancer biology across human, mouse, and rat samples.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FGFR2 Rabbit mAb
Remarks/Alias	BEK; JWS; BBDS; CEK3; CFD1; ECT1; KGFR; TK14; TK25; BFR-1; CD332; K-SAM; R2
Species	Human
GeneID (Human)	2263
GeneID	2263
Immunogen	Synthetic Peptide: A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGFR2 (P21802).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21802
Protein Weight	92kDa
Shipping	Ice bag