


Anti-Fgfr2 Polyclonal Antibody For Wb, If/Icc, Elisa - P21802

Item Number: CM0013129



Introduction

Rabbit polyclonal antibody against human FGFR2, validated for WB, IF/ICC, ELISA. Reacts with human, mouse, rat. Ideal for studying cell proliferation, differentiation, and developmental signaling pathways.

[Learn More](#)

Parameter	Value
Product Name	FGFR2 Rabbit pAb
Remarks/Alias	BEK; JWS; BBDS; CEK3; CFD1; ECT1; KGFR; TK14; TK25; BFR-1; CD332; K-SAM; FGFR2
Species	Human
GeneID (Human)	2263
GeneID	2263
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P21802
Protein Weight	92kDa
Shipping	Ice bag