

Anti-Fgfr3 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P22607

Item Number: CM0017553



Introduction

Rabbit polyclonal antibody targeting human FGFR3, specified for WB, IHC-P, IF/ICC, and ELISA. Cross-reacts with mouse and rat. An essential tool for skeletal development and FGFR3-related cancer research.

[Learn More](#)

Parameter	Value
Product Name	FGFR3 Rabbit pAb
Remarks/Alias	ACH; CEK2; JTK4; CD333; H5FGFR3EX; FGFR3
Species	Human
Gene ID (Human)	2261
Gene ID	2261
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 700-806 of human FGFR3 (NP_000133.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22607
Protein Weight	88kDa
Shipping	Ice bag