

Anti-Fgfr3 Rabbit Monoclonal Antibody For Ihc-P And Elisa - P22607

Item Number: CM0003384



Introduction

FGFR3 rabbit monoclonal antibody (mAb) validated for IHC-P and ELISA applications. Specifically recognizes human FGFR3 (P22607). Ideal for investigating skeletal development, cancer signaling, and receptor tyrosine kinase pathways.

[Learn More](#)

Parameter	Value
Product Name	FGFR3 Rabbit mAb
Remarks/Alias	ACH; CEK2; JTK4; CD333; HSGFR3EX; FGFR3
Species	Human
Gene ID (Human)	2261
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FGFR3 (P22607)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	P22607
Protein Weight	88kDa
Shipping	Ice bag