


Anti-Fgfr4 Monoclonal Antibody For Wb, Elisa - P22455

Item Number: CM0004201



Introduction

Recombinant rabbit monoclonal antibody against human FGFR4 (P22455). Suitable for Western blot and ELISA, with cross-reactivity to mouse and rat. Ideal for research in cell signaling, lipid metabolism, and oncology.

[Learn More](#)

Parameter	Value
Product Name	FGFR4 Rabbit mAb
Remarks/Alias	TKF, JTK2, CD334, FGFR4
Species	Human
GeneID (human)	2264
GeneID	2264
Immunogen	Synthetic peptide (aa 700-802 of human FGFR4, P22455)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22455
Protein Weight	88 kDa
Shipping	Ice bag