


Anti-Fgfr4 Polyclonal Antibody For Wb, Elisa - P22455

Item Number: CM0028835



Introduction

High-quality rabbit polyclonal antibody against human FGFR4 for Western blotting and ELISA. Recognizes human, mouse, and rat FGFR4 protein. Ideal for investigating FGFR4 signaling in cell proliferation, metabolism, and cancer.

[Learn More](#)

Parameter	Value
Product Name	FGFR4 Rabbit pAb
Remarks/Alias	TKF; JTK2; CD334; FGFR4
Species	Human
Gene ID (Human)	2264
Gene ID	2264
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 22-150 of human FGFR4 (NP_998812.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22455
Protein Weight	88kDa
Shipping	Ice bag