

Anti-Fgg Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - P02679

Item Number: CM0009763



Introduction

Recombinant rabbit monoclonal antibody targeting human Fibrinogen gamma chain (FGG). Validated for WB, IF/ICC, IHC-P, and ELISA; detects ~49/52 kDa bands. Ideal for thrombosis and hemostasis research.

[Learn More](#)

Parameter	Value
Product Name	FGG Rabbit mAb
Remarks/Alias	FGG
Species	Human
Gene ID (Human)	2266
Gene ID	2266
Immunogen	Recombinant protein
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	P02679
Protein Weight	49 kDa/52 kDa
Shipping	Ice bag