

Anti-Fgg Rabbit Polyclonal Antibody For Wb And Elisa - P02679

Item Number: CM0026190



Introduction

Rabbit polyclonal antibody against human FGG (Fibrinogen gamma chain), suitable for WB and ELISA. Cross-reacts with mouse. Ideal for hemostasis and blood coagulation research.

[Learn More](#)

Parameter	Value
Product Name	FGG Rabbit pAb
Remarks/Alias	FGG
Species	Human
GeneID (human)	2266
GeneID	2266
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-330 of human FGG (NP_000500.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02679
Protein Weight	52kDa
Shipping	Ice bag