

Anti-Fgl2 Polyclonal Antibody For Wb And Elisa - Q14314

Item Number: CM0015726

Introduction

High-quality rabbit polyclonal antibody targeting human Fibroblast Growth Factor 2 (FGF2). Validated for Western blot and ELISA. Suitable for exploring lymphocyte functions at mucosal sites.

[Learn More](#)

Parameter	Value
Product Name	FGL2 Rabbit pAb
Remarks/Alias	FGL2; T49; pT49
Species	Human
GeneID (human)	10875
GeneID	10875
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 24-140 of human FGL2 (NP_006673.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q14314
Protein Weight	50kDa
Shipping	Ice bag