

Anti-Fh (Fumarate Hydratase) Polyclonal Antibody For Wb, Elisa - P07954

Item Number: CM0022516



Introduction

Rabbit polyclonal antibody against human FH, validated for WB and ELISA. Recognizes human, mouse, rat FH. Suitable for studies of TCA cycle, DNA repair, and tumor suppression.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] FH Rabbit pAb
Remarks/Alias	MCL; FMRD; HsFH; LRCC; HLRCC; MCUL1; FH
Species	Human
Gene ID (Human)	2271
Gene ID	2271
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 50-150 of human FH (NP_000134.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P07954
Protein Weight	55kDa
Shipping	Ice bag