

## Anti-Fh (Fumarate Hydratase) Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P07954

Item Number: CM0027157



### Introduction

Rabbit polyclonal antibody against human FH (Fumarate Hydratase), validated for WB, IHC-P, IF/ICC, ELISA, with cross-reactivity to mouse and rat. Suitable for TCA cycle and DNA repair studies.

[Learn More](#)

| Parameter        | Value                                                                                                           |
|------------------|-----------------------------------------------------------------------------------------------------------------|
| Product Name     | FH Rabbit pAb                                                                                                   |
| Remarks/Alias    | MCL; FMRD; HsFH; LRCC; HLRCC; MCUL1; FH                                                                         |
| Species          | Human                                                                                                           |
| Gene ID          | 2271                                                                                                            |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 45-510 of human FH (NP_000134.2). |
| Source           | Rabbit                                                                                                          |
| Category         | Polyclonal Antibodies                                                                                           |
| Application      | WB, IHC-P, IF/ICC, ELISA                                                                                        |
| Cross Reactivity | Human, Mouse, Rat                                                                                               |
| SWISS            | P07954                                                                                                          |
| Protein Weight   | 55kDa                                                                                                           |
| Shipping         | Ice bag                                                                                                         |