

Anti-Fhit Monoclonal Antibody For Wb, Elisa - P49789

Item Number: CM0005973



Introduction

Rabbit monoclonal antibody targeting human FHIT (FRA3B/AP3Aase). Validated for WB and ELISA. Cross-reacts with human, mouse, and rat. 17kDa protein. For apoptosis and tumor suppressor research.

[Learn More](#)

Parameter	Value
Product Name	FHIT Rabbit mAb
Remarks/Alias	FRA3B; AP3Aase; FHIT
Species	Human
GeneID (Human)	2272
GeneID	2272
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human FHIT (P49789).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P49789
Protein Weight	17kDa
Shipping	Ice bag