

Anti-Fhit Polyclonal Antibody For Wb, Elisa - P49789

Item Number: CM0015052



Introduction

Rabbit polyclonal antibody against human FHIT (FRA3B, AP3Aase), validated for Western blot and ELISA. Suitable for studying tumor suppression, apoptosis, and nucleotide metabolism in human and mouse samples.

[Learn More](#)

Parameter	Value
Product Name	FHIT Rabbit pAb
Remarks/Alias	FRA3B, AP3Aase, FHIT
Species	Human
GeneID (Human)	2272
GeneID	2272
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 25-147 of human FHIT (NP_001159715.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P49789
Protein Weight	17kDa
Shipping	Ice bag