

Anti-Fhl1 Monoclonal Antibody For Ihc-P, Elisa - Q13642

Item Number: CM0011929



Introduction

High-quality rabbit monoclonal antibody against human FHL1, validated for IHC-P and ELISA applications. Suitable for muscle development and disease research. Specific to human FHL1 (Q13642).

[Learn More](#)

Parameter	Value
Product Name	FHL1 Rabbit mAb
Remarks/Alias	KYOT, SLIM, FCMSU, FHL-1, FHL1A, FHL1B, FLH1A, SLIM1, XMPMA, RBMX1A, RBMX1B, SLIM-1, SLIMMER, FHL1
Species	Human
GeneID	2273
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human FHL1 (Q13642)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q13642
Protein Weight	36kDa
Shipping	Ice bag