


Anti-Fhl2 Monoclonal Antibody For Wb, Elisa - Q14192

Item Number: CM0004702



Introduction

Rabbit monoclonal antibody targeting human FHL2, validated for WB and ELISA. Cross-reacts with mouse and rat. Ideal for studying transcriptional regulation, apoptosis, and muscle biology. Q14192.

[Learn More](#)

Parameter	Value
Product Name	FHL2 Rabbit mAb
Remarks/Alias	DRAL; AAG11; FHL-2; SLIM3; SLIM-3; FHL2
Species	Human
GeneID (Human)	2274
GeneID	2274
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 50-260 of human FHL2 (Q14192).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q14192
Protein Weight	32kDa
Shipping	Ice bag