


Anti-Fhl3 Polyclonal Antibody For Wb, If/Icc, Elisa - Q13643

Item Number: CM0030756



Introduction

Rabbit anti-FHL3 polyclonal antibody validated for WB, IF/ICC, and ELISA. Detects a 31kDa protein in human and mouse skeletal muscle, supporting transcription coactivator research.

[Learn More](#)

Parameter	Value
Product Name	FHL3 Rabbit pAb
Remarks/Alias	SLIM2; FHL3
Species	Human
GeneID (Human)	2275
GeneID	2275
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13643
Protein Weight	31kDa
Shipping	Ice bag