

Anti-Fhl5 Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - Q5Td97

Item Number: CM0008455



Introduction

High-affinity rabbit monoclonal antibody against human FHL5 (ACT), validated for WB, IF/ICC, IP, and ELISA. Cross-reacts with mouse. Suitable for spermatogenesis and testis-specific transcriptional regulation research.

[Learn More](#)

Parameter	Value
Product Name	FHL5 Rabbit mAb
Remarks/Alias	ACT; FHL-5; dj393D12.2; 1700027G07Rik; FHL5
Species	Human
GeneID	9457
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 162-284 of human FHL5 (NP_065228.4)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5TD97
Protein Weight	33 kDa
Shipping	Ice bag