

Anti-FHL5 Rabbit Polyclonal Antibody For Wb, Elisa - Q5Td97

Item Number: CM0017070



Introduction

High-quality rabbit polyclonal antibody targeting human FHL5 (ACT), suitable for Western blot and ELISA applications. Detects human and mouse FHL5. Ideal for spermatogenesis and CREM transcriptional activity studies.

[Learn More](#)

Parameter	Value
Product Name	FHL5 Rabbit pAb
Remarks/Alias	ACT; FHL-5; dj393D12.2; 1700027G07Rik; FHL5
Species	Human
GeneID	9457
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-284 of human FHL5 (NP_065228.4).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q5TD97
Protein Weight	33kDa
Shipping	Ice bag