

Anti-Fibp Rabbit Polyclonal Antibody For Wb, Ihc-P - O43427

Item Number: CM0025974

Introduction

Rabbit polyclonal antibody against human FIBP (Acidic FGF intracellular-binding protein) for WB, IHC-P, IP, ELISA. Validated in human, mouse, rat. Involved in FGF1 signaling and embryonic laterality.

[Learn More](#)

Parameter	Value
Product Name	FIBP Rabbit pAb
Remarks/Alias	FGFIBP; FIBP-1; TROFAS; FIBP
Species	Human
GeneID (Human)	9158
GeneID	9158
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-364 of human FIBP (NP_004205.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O43427
Protein Weight	42kDa
Shipping	Ice bag