


Anti-Fibrillarin Rabbit Monoclonal Antibody For Wb - P22087

Item Number: CM0005659



Introduction

Rabbit monoclonal antibody against human Fibrillarin (FBL, 34 kDa), validated for WB, IP, IF/ICC, IHC-P, ELISA. Reacts with human, mouse, rat. Ideal for nucleolar protein and rRNA methylation research.

[Learn More](#)

Parameter	Value
Product Name	Fibrillarin Rabbit mAb
Remarks/Alias	FIB; FLRN; Nop1; RNU3IP1
Species	Human
GeneID (Human)	2091
GeneID	2091
Immunogen	Synthetic peptide: A synthetic peptide corresponding to a sequence within amino acids 111-210 of human Fibrillarin/U3 RNP (NP_001427.2).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IP, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22087
Protein Weight	34kDa
Shipping	Ice bag