

Anti-Fibrillarin/U3 Rnp Rabbit Monoclonal Antibody For Wb, If/Icc, Ip, Elisa - P22087

Item Number: CM0003045



Introduction

Rabbit monoclonal antibody against human Fibrillarin/U3 RNP (P22087), validated for WB, IF/ICC, IP, ELISA. Detects 34 kDa protein; cross-reacts with mouse and rat. Suitable for nucleolar biology studies.

[Learn More](#)

Parameter	Value
Product Name	Fibrillarin/U3 RNP Rabbit mAb
Remarks/Alias	FIB; FLRN; Nop1; RNU3IP1; Fibrillarin/U3 RNP
Species	Human
Gene ID (Human)	2091
Gene ID	2091
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 222-321 of human Fibrillarin/U3 RNP (P22087).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P22087
Protein Weight	34kDa
Shipping	Ice bag